

# Ocean Optics USB4000 with *usb4java*

Bernard Panneton  
December 2019

## Introduction

The `usb4java` library<sup>1</sup> is used to communicate with a USB peripheral. It used to interact with an Ocean Optics Inc. USB4000 spectrometers.

This library was used to build a collection of R functions stored in ***playWith\_\_usbjava.R*** file. The commands required to communicate with the spectrometer are detailed in the device technical manual<sup>2</sup>.

## Functions in `playWith__usb4java.R`

### `init__usb()`

Initialize usb device:

1. load required library
2. init JVM
3. Set Java class paths
4. Define some objects

RETURN: a list of 4 components:

1. context: a Java object of class `org.usb4java.Context`
2. dlist: a Java object of class `org.usb4java.DeviceList`
3. libusb: a Java object of class `org.usb4java.LibUsb`
4. bufutils: a Java object of class `org.usb4java.BufferUtils`

### `find__usb <- function(product,vendor,usbObjects, silent=TRUE)`

Given a vendor and a product ID number, find the corresponding USB device.

INPUTS:

- product: product ID number
- vendor: vendor ID number
- usbObjects: list returned by `init__usb`
- silent: when TRUE, no output at console.

OUTPUTS:

- a list with 2 components
  1. `usbDevice`: the device as obtained with `dlist$get(as.integer(dev_no))` where `dlist` was defined in `init__usb`

---

<sup>1</sup><http://usb4java.org/>

<sup>2</sup>USB4000-OEM-Data-Sheet.pdf stored in the **Doc** of the RStudio project **OceanOptics\_\_with\_\_usb4java\_in\_R**.

2. usbDescription: obtained by `libusb$DeviceDescriptor(usbDevice, usbDescription)`

### **get\_OO\_name\_n\_serial(usbObjects, usbDevice)**

Function to get device name and version, device serial number and company names.

INPUTS:

- usbObjects: the list returned by `init_usb()`
- usbDevice: the list returned by `find_usb()`

OUTPUTS:

- list with 3 components:
  1. name: name of the USB device
  2. version: name of device with version number
- 3. serialno: serial number

### **getWavelengths(usbObjects, usbDevice){**

To get the wavelength vector by reading the wavelength calibration coefficients

INPUTS:

- usbObjects: the list returned by `init_usb()`
- usbDevice: the list returned by `find_usb()`

OUTPUTS:

- a vector of wavelengths

### **jbyte\_2\_uint(x)**

Takes a vector of Java bytes and interprets and 0:255

INPUTS:

- x: vector of Java bytes as seen in R

OUTPUTS:

- a vector of value in the range 0:255, same length as input.

### **queryStatus <- function(usbObjects, usbDevice)**

Query the USB device status.

INPUTS:

- usbObjects: the list returned by `init_usb()`
- usbDevice: the list returned by `find_usb()`

OUTPUTS:

- a list of 5 elements:
  1. nb\_pix: number of pixels in spectrum
  2. int\_time: current integration time
  3. pack\_in\_spectra: number of data packets per spectrum
  4. pack\_count:
  5. usb\_speed: speed of USB transfer (“full” or “high”)

### **revShort\_2\_numeric <- function(x)**

Function to reverse byte order and generate numeric vector.

INPUTS:

- x: a vector of bytes as returned by a Java byte buffer interpreted as an array using `.jevalArray()`

### **OUTPUTS:**

- a numeric vector where each value is from a pair of bytes grouped in reverse order to build a word that is converted to numeric.

### **setIntegrationTime(temps,usbObjects, usbDevice)**

Function to set spectrometer integration time.

INPUTS:

- temps: integration time in msec.
- usbObjects: the list returned by `init_usb()`
- usbDevice: the list returned by `find_usb()`

OUTPUTS:

- none

### **getSpectrum(pack\_in\_spectra=15, usbObjects, usbDevice)**

Function to retrieve a spectrum.

INPUTS:

- pack\_in\_spectra: number of data packets per spectrum
- usbObjects: the list returned by `init_usb()`
- usbDevice: the list returned by `find_usb()`

OUTPUTS:

- a spectrum as a numeric vector.

### **get\_N\_Spectrum(pack\_in\_spectra=15, nspectra=2, usbObjects, usbDevice)**

Function to retrieve a spectrum made as an average over a number of spectra.

INPUTS:

- pack\_in\_spectra: number of data packets per spectrum
- nspectra: number of spectrum to average over.
- usbObjects: the list returned by `init_usb()`
- usbDevice: the list returned by `find_usb()`

OUTPUTS:

- a spectrum as a numeric vector.

## free\_Device(usbObjects)

Function to free device. Required to cleanly end.

INPUTS:

- – usbObjects: the list returned by init\_usb()

## Example

Connect a USB4000 to a USB port and copy the following lines to the R console.

```
#Testing functions
source("R/playWith_usb4java.R")
usbObjects <- init_usb()
product=0x1022 vendor=0x2457 usbDevice <- find_usb(product,vendor,usbObjects,TRUE)
name_serial <- get_OO_name_n_serial(usbObjects, usbDevice$usbDevice) lapply(name_serial, print)
wv <- getWavelengths(usbObjects, usbDevice$usbDevice)
statut <- queryStatus(usbObjects, usbDevice$usbDevice)
setIntegrationTime(70,usbObjects,usbDevice$usbDevice)(statut <- queryStatus(usbObjects,usbDevice$usbDevice))
dum <- getSpectrum(pack_in_spectra=15, usbObjects, usbDevice$usbDevice)dum <- getSpectrum(pack_in_spectra =
15,usbObjects,usbDevice$usbDevice)
plot(wv, dum[22:3669],type="l",col="red",lwd=2, ylim=c(0,7000))
sp <- get_N_Spectrum(pack_in_spectra=15, nspectra=20, usbObjects, usbDevice$usbDevice)plot(wv, sp[22 :
3669], type = "l", col = "red", lwd = 2, ylim = c(0, 7000), main = paste0(name_serial$name, " - Serial number:
", name_serial$serialno), xlab = "Wavelength [nm]", ylab = "Intensity [A.U.]")
free_Device(usbObjects)
```

## Final words

The code in the R project works with a USB4000. Adaptation to other spectrometer should be fairly straightforward. One needs to check the product ID and the exact syntax of the various commands and some parameters and the number of pixels in a spectrum, the number of packets transmitted over the USB...

```
rm(list=ls())
```